

HEEL PILLOW - 100% cotton



Product Code	28556
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2386673
CND	V0811
EAN/UPC	8023279285567
PARAF	975385978

Description

HEEL CUSHION

Size: 32 x 12 cm

Products to prevent low decubitus risk, made of silicone fibre, coated with 100% preshrunk cotton cover. Main property of silicone fibre is to maintain original shape also after strong pressure. This allows to maintain comfort and proper weight distribution in the time. Breathable cotton cover provides maximum cool comfort. Washable at 30-40°C.

All cushions and pillow are supplied with multilingual instructions: GB, FR, IT, ES, PT, DE

Made in Italy.