



Product Code	28565
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2398246
CND	V080701
EAN/UPC	8023279285659
PARAF	975386006

## Description

Special silent pump for an excellent pressure regulation by a knob on the front panel.

It can be hung to any kind of hospital bed frame by means of two adjustable hooks.

For use with bubble or alternating cell mattresses.

## Technical Specifications

Size: 25x12x10 cm  
 Operating voltage: 110-220 V, 50-60 Hz  
 Working pressure: 40-100 mmHg  
 Air output: 6-7 liter/min  
 Cycle time: 6 min  
 Case: TUV-CE IEC/EN60601-1; 60601-1-2