


INTERCHANGEABLE CELL AIR MATTRESS + COMPRESSOR (28567+28568)



Product Code	28569
Unit of sale	1 kit
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2398257
CND	V080701
EAN/UPC	8023279285697
PARAF	975386032

Description

KIT INTERCHANGEABLE CELL-MATTRESS + PUMP SYSTEM

• INTERCHANGEABLE CELL-MATTRESS

Recommended for the prevention and treatment of pressure sore stage I and II in medium-term and domiciliary therapies. The 20 9.5 cm-high Nylon+PVC cells provide superior comfort. Cells (20+1 spare) are fixed by means of pressure stud and can easily be replaced in case of breakage. Manual and box in 10 languages: GB, FR, IT, ES, DE, PT, PL, HU, GR, Arabic.

• PUMP SYSTEM

Special silent pump for an excellent pressure regulation by a knob on the front panel. It can be hung to any kind of hospital bed frame by means of two adjustable hooks. For use with bubble or alternating cell mattresses.

Technical Specifications

Bubble mattress:

Total size: 200 x 86 x 9.5 cm
Cell height: 9.5 cm (4")
Material: Nylon+PVC
Thickness: 0.32 mm
Suitable for persons up to: 140 kg

Pump system:

Size: 25x12xh 10 cm
Operating voltage: 110-220 V, 50-60 Hz
Working pressure: 40-100 mmHg
Air output: 6-7 liter/min
Cycle time: 12 min
Case: TUV-CE IEC/EN60601-1; 60601-1-2