

KIT 5010E ALTERNATING PRESSURE MATTRESS + PUMP



Product Code	28575
Unit of sale	1 kit
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2472448
CND	V080701
EAN/UPC	8023279285758
PARAF	975386057

Description

KIT 5010E ALTERNATING PRESSURE MATTRESS + PUMP

MODERATE TO INTENSIVE CARE STAGE I, II, III

Valuable overlay system designed for prevention and treatment of patients at moderate to intensive risk of pressure ulcers. Pressure is periodically redistributed by inflating and deflating cells in order to mimic natural sleep movement patterns and avoid long term pressurization of tissue.

Manual in: GB, FR, IT, ES, PT, DE, PL, RO, SE, GR.

Ventilated low-air-loss cells wicks moisture away

Professional use

Nylon PU cover with button biocompatibility, vapor permeable & water resistant. Flame Retardant.

ALTERNATING FUNCTION

The pump inflates and deflates the cells with cyclicity, to change the contact position between the patient and the bed, to promote the blood circulation, to prevent and recover the bedsore, reducing the patient's pain and the nursing labor.

Technical Specifications

PUMP:

Size: 30x13.5x11 cm
Cycle time: 12 min
Operating voltage: 220-240 V, 50 Hz
Working pressure: 55-120 mmHg
Case material: Flame retardant ABS
Air output: 7-8 liter/min
TUV-CE IEC/EN60601-1; 60601-1-2

MATTRESS:

Total size: 198x86x12.8 cm
Alternating 1 in 2 (group A and B)
Cell height: 12.8 cm (5")
Cover material: Nylon PU with button
Thickness: 0.4 mm
Number of cells: 17+1 spare
Cells material: nylon PVC
Load: 150 kg