

ICE BAG diameter 28 cm - large cap 6 cm



Product Code	28605
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2831166
CND	M9001
EAN/UPC	8023279286052
PARAF	925813851
GTIN	08023279286052
GMDN	12042

Description

ICE BAG - Ø 28 cm - large top cap - light blue

Cap Ø: external 8.78 cm, internal 6 cm

Single packed in carton box, with multilingual manual (GB, FR, IT, PL, ES, PT, GR, Arabic).

Made of polyester fabric + PVC.

Neck of bag: PP