

HEATING UNDERBLANKET - single



Product Code	28660
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279286601
PARAF	950453946

Description

HEATING UNDERBLANKET - single

Electric underblankets made of non-woven fabric (polyester) with stamp welding manufacturing. They warm and relax the entire body, night and day.
3 temperature settings with level indicator.
Detachable hand controller allows to wash the blanket.
Overheating protection.
Supplied with multilingual box and manual GB, FR, IT, DE, NL, PT, GR, SE, PL, RO, FI, CZ.

Technical Specifications

Size (cm): 150x80
Power consumption: 60 W
Voltage: 220-240 V