

HEATING UNDERBLANKET - double



Product Code	28661
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279286618
PARAF	950454076

Description

HEATING UNDERBLANKET - double

Electric underblankets made of non-woven fabric (polyester) with stamp welding manufacturing. They warm and relax the entire body, night and day.
3 temperature settings with level indicator.
Double underblanket has 2 controllers.
Detachable hand controller allows to wash the blanket.
Overheating protection.
Supplied with multilingual box and manual GB, FR, IT, DE, NL, PT, GR, SE, PL, RO, FI, CZ.

Technical Specifications

Size (cm): 160x140
Power consumption: 2x60 W
Voltage: 220-240 V