

# HEATING OVERBLANKET 160x120 cm



Product Code	28664
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279286649
PARAF	950454090

## Description

Electric over blanket made of soft micro plush provides warmth and comfort night and day.  
3 temperature setting with LED indicator (0, 1, 2, 3).  
Automatic switch off after 3 hours.  
Detachable hand controller allows to wash the blanket.  
Fast heating and overheating protection.  
Supplied with multilingual box and manual GB, FR, IT, DE, NL, PT, GR, SE, PL, RO, FI, CZ.

## Technical Specifications

Size (cm): 160x120  
Power consumption: 120 W  
Voltage: 220-240 V