

# PEDIATRIC URINE BAG 100 ml - sterile



Product Code	28685
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	Yes

RDM (NSIS)	1426154
CND	A06030301
EAN/UPC	8023279286854
PARAF	976329829

## Description

Paediatric urine bag, individually wrapped in polybag.  
They can be used for boy and girl.  
Graduated 10 to 100 ml, sterile.

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, SE, FI, PL, CZ, HU, HR, RO, GR.