

URINE BAG 2000 cc - tube 90 cm - with non-return/exhaust valve



Product Code	28691
Unit of sale	box of 30 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2205324
CND	A06030399
EAN/UPC	8023279286915
PARAF	970523751

Description

URINE BAGS WITH NON RETURN VALVE

Urine bag 2000 cc - tube 90 cm - polybag

Minimum order: box of 30*

*Polybag of 30 bags with tips (Italian reimbursement label) and 1 hanger

Instructions in GB, FR, IT, ES, DE, PT, GR, PL. Non sterile.

Tube Ø: internal 5.05 mm, external 6.35 mm. Latex and DEHP free.

Made in Spain.