



Product Code	28720
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2339132
CND	G020301
EAN/UPC	8023279287202
PARAF	973923004
GMDN	11582

Description

ENEMA IRRIGATOR - 2 l capacity

Supplied in cardboard box with 2 cannulas (vaginal, anal), tap with opening lever, tube and carrying bag.

Multilanguage box: GB, FR, IT, ES.