

# Z STRAP 3 in 1 FIXING SYSTEM for urine bag, tube, catheter



Product Code	28727
Unit of sale	box of 2 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2491623
CND	T0499
EAN/UPC	7290018953600
PARAF	950080349

## Description

Comfortable and discreet all in one universal strap for fixing a urine bag, a urine bag tube or a catheter.

Special designed Z-strap is made of top-quality medical garment preventing irritations; silicone beading on internal side keeps the strap in place preventing any inconvenience.

- one size fits all (70 cm long), can be cut to desired size
- wide band (5 cm width), distributes pressure and weight better, especially for diabetics
- easy application, for nurse and user, fast wear and removal
- skin friendly, no use of adhesive / no need to shave hair
- fixing straps for drainage tubes to prevent slipping
- prevents injuries caused by sudden movement
- fits all size of urine bags / leg bags
- reusable can be machine washed, does not shrink
- prevents pressure wounds, wide strap less pressure per cm<sup>2</sup>