

MEN'S EXTERNAL CATHETER with tray - silicone - diam 31 mm



Product Code	28737
Unit of sale	box of 30 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2671376
CND	U070201
EAN/UPC	6946680199020
PARAF	950481376

Description

MALE EXTERNAL SILICONE CATHETERS

- Male ext. catheter with tray Ø 31 mm
- Size: M/L

Male external catheter with tray made of 100% ultrathin silicone, that eliminates risk of latex related skin irritations. This material has the highest water vapor permeability. Soft, breathable and odour free silicone provides greater patient's comfort and healthier skin with less risk of maceration. Adhesive band is part of the catheter sheath: it facilitates the application and offers patient's safety and extended wearing period. Integrated tray ensures that adhesive ends do not stick together. Its transparency allows easy skin inspection.

Security

Kink-resistant silicone funnel helps to ensure continuous urine flow and secured connection to all major types of urine collectors.

Ease of use

One-piece system easy to apply and to remove, secured to the skin with an acrylic adhesive.

A perfect fit for most patients

It has a standard-length sheath with strong adhesive and is available in five different sizes.

Individually packed in a plastic bag in a box of 30.