

PREMIUM INFANT CUFF + 1 TUBE TPU BLADDER



Product Code	32886
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2512518
CND	C99
EAN/UPC	8023279328868
GTIN	08023279328868
GMDN	34978

Description

PREMIUM CUFFS dark grey - nylon - TPU bladder - velcro fastening
Made of high quality nylon 210D.
Durability test up to 10,000 cycles.

- One tube cuff for Sirio, London, Boston, Roma 1T, Tokyo, Madrid, Ri-san, Precisa, Minimus...
- Model: Infant
- Size cm: 7.5x26
- Arm circumference: 13 - 19 cm