

PREMIUM INFANT CUFF + 2 TUBE TPU BLADDER



Product Code	32887
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2512518
CND	C99
EAN/UPC	8023279328875
PARAF	987471404
GTIN	08023279328875
GMDN	34978

Description

PREMIUM CUFFS dark grey - nylon - TPU bladder - velcro fastening
Made of high quality nylon 210D.
Durability test up to 10,000 cycles.

Two tubes cuff for Minor, Yton line, Roma 2T, Dallas, Big Ben, Empire, Dayton...

- Model: Infant
- Size cm: 7.5x26
- Arm circumference: 13 - 19 cm