

TRUBARB ABSORB. SUTURE gauge 3/0 circle 3/8 needle
19mm - 45cm - violet



Product Code	22896
Unit of sale	box of 6 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2380161
CND	H0101010101
EAN/UPC	8023279228960

Description

TRUBARB - absorbable surgical needed-suture (synthetic) knotless tissue-closure device

Equivalent Ethicon brand: STRATAFIX

- Gauge: 3/0
- Metric: 2
- Colour: Violet
- Circle: 3/8
- Needle profile: Reverse cutting
- Needle length: 19 mm
- Suture length: 45 cm

- Technical details
- monofilament polydioxanone
 - secure and effective triangular end-stopper
 - single angle cut uni-directional helical barb design
 - world class 300 series QNPL needle
 - sterilized by ethylene oxide
 - shelf life of 3 years

- Applications
- gynaecological and urology procedures
 - general surgery procedures
 - oncology procedures
 - plastic and cosmetic surgeries

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.