



Product Code	28972
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2535067
CND	Y181899
EAN/UPC	8023279289725
PARAF	925875320
GMDN	35584

## Description

### SAFETY GRAB BAR - 443 mm

Innovative line of GRAB BAR, can be quickly installed everywhere without need of drilling and screwing.  
 Patented suction cups (Ø 93 mm).  
 Holds up to 60 kg vertical and 45 lateral weight (tested by SGS).

Ideal solution to support elder people during daily bathroom use.  
 Top quality product made of heavy duty ABS, with patented "green/red" suction indicator for better safety.

Multilingual box and manual: GB, FR, IT, ES, PT, DE, GR, Arabic.