

PREGNANCY TEST - self test - midstream (large wipe) - 1 test



Product Code	29103
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

UK-REP	No
CH-REP	No

RDM (NSIS)	1702289
CND	W0102160302
EAN/UPC	8023279291032
PARAF	977812890

Description

PREGNANCY TEST

- SELF TEST
- Type: midstream (large wipe)
- Package: box of 1 pc

GIMA offers a complete line of pregnancy tests, both for OTC market (self test) and for professional use in hospital or doctor's cabinets. All products are packed in cardboard boxes with multilingual user manual (GB, IT, ES, FR, GR, PT).

Main features of pregnancy tests:

- each test is individually packed in sealed pouch
- very quick: results in 3 minutes
- early detection: it can be used since the first day you miss your period
- accuracy: 99%
- sensitivity: 10 mIU/ml
- Clinically tested

SHELF LIFE 36 MONTHS