



Product Code	32946
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2556483
CND	Z1203020501
EAN/UPC	401567211141
PARAF	976386122

Description

OMRON RS4 HEM-6181-E - wrist

LED positioning sensor and cuff wrap guide indicate if the cuff is positioned correctly.

Multilingual box and manual: GB, FR, IT, ES, DE, NL, RU, TR, Arabic.

Technical Specifications

Measuring method: oscillometric

Measuring range:

Pressure: 0 mmHg to 299 mmHg $\pm$ 3 mmHg

Pulse: 40 to 180 min $\pm$ 5% of display reading

Memory (with date and time): 60

Inflation: automatic

Deflation: automatic

Rating: DC 3 V 3.0 W

Power source: 2 "AAA" batteries 1.5 V

Battery life (No. of measurements): 300

Protection against electric shock: internally powered ME equipment (when using batteries)

Weight (console + cuff) without batteries: 86 g

Size (mm): 93 $\times$ 62 $\times$ h 20

Cuff materials: nylon and polyester

Measurable circumference (cm): 13.5-21.5 (wrist)