

GIMASONIC POCKET DOPPLER WITH COLOUR DISPLAY without probe



Product Code	29480
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2681945
CND	Z12080103
EAN/UPC	8023279294804
UMDNS	11696

Description

COLOUR SCREEN POCKET DOPPLER

GIMASONIC Colour screen doppler - without probe

- built-in speaker
 - TFT-LCD, real-time display of FHR
 - three working modes: real-time FHR, average FHR and manual calculation
 - display of FHR, bar graph and heartbeat waveform, red alarm for abnormal FHR
 - low power consumption, continuously work for more than 8 hours with two "AA" batteries
 - indication for probe-off and low battery
 - range of probes for foetal or vascular use: 2, 3, 8 MHz probe
- Suitable for use after the 12th week of pregnancy
Multilingual box and manual: GB, FR, IT, ES, PT, DE.

Technical Specifications

- Ultrasound working frequency: 2 MHz $\pm$ 10%
- Fetal heart rate:
 - range: 50-240 bpm
 - resolution: 1 bpm
 - accuracy: $\pm$ 2 bpm
- Alarm:
 - high limit: 160 bpm
 - low limit: 120 bpm
- Ultrasound intensity: < 5 mW/cm²
- Auto power off after 1 minute no signal
- Dimension: 135x92x29 mm
- Weight: 245 g
- Power supply: 2 AA Batteries LR6 - 1.5 V DC