



Product Code	29490
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

CH-REP	No
--------	----

RDM (NSIS)	2958300
CND	Z12080103
EAN/UPC	6944413800120
UMDNS	11696
GMDN	34040

Description

SONOTRAX LITE DOPPLER - without probe - with bag

Technical Specifications

Size: 32x85x138 mm
 Weight: 290 g (including battery)
 Heart rate range: 50-210 bpm
 Resolution: 1 bpm
 Accuracy: ± 3 bpm
 LCD display: 45x25 mm (only code 29491-93)
 Auto average and manual mode
 Recording length: 480 s / 240 s
 Auto power-off to save power
 Manual: GB, IT - on request: DE, ES, FR, PT.