

SONOTRAX BASIC POCKET DOPPLER WITH DISPLAY without probe



| | |
|---------------|----------------|
| Product Code | 29491 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|------------|---------------|
| CH-REP | No |
| RDM (NSIS) | 2755473 |
| CND | Z12059005 |
| EAN/UPC | 6944413802384 |
| UMDNS | 11696 |
| GMDN | 34040 |

Description

SONOTRAX BASIC DOPPLER - with display - without probe - with carrying bag

Technical Specifications

Size: 32x85x138 mm
Weight: 290 g (including battery)
Heart rate range: 50-210 bpm
Resolution: 1 bpm
Accuracy: ± 3 bpm
LCD display: 45x25 mm (only code 29491-93)
Auto average and manual mode
Recording length: 480 s / 240 s
Auto power-off to save power
Manual: GB, IT - on request: DE, ES, FR, PT.