

# SONOTRAX PRO POCKET DOPPLER WITH DISPLAY

## without probe



Product Code	29493
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

CH-REP	No
RDM (NSIS)	2640520
CND	Z12059005
EAN/UPC	6944413800175
UMDNS	11696
GMDN	34040

### Description

SONOTRAX PRO II DOPPLER with backlighted display built-in audio recorder and rechargeable batteries with battery charger - without probe - with carrying bag

### Technical Specifications

Size: 32x85x138 mm  
Weight: 290 g (including battery)  
Heart rate range: 50-210 bpm  
Resolution: 1 bpm  
Accuracy:  $\pm 3$  bpm  
LCD display: 45x25 mm (only code 29491-93)  
Auto average and manual mode  
Recording length: 480 s / 240 s  
Auto power-off to save power  
Manual: GB, IT - on request: DE, ES, FR, PT.