

# SONOTRAX 4 MHz PROBE



|               |                |
|---------------|----------------|
| Product Code  | 29498          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|        |    |
|--------|----|
| CH-REP | No |
|--------|----|

|            |              |
|------------|--------------|
| RDM (NSIS) | 213512       |
| CND        | Z1208010480  |
| EAN/UPC    | 694441380552 |
| UMDNS      | 11696        |
| GMDN       | 34040        |

## Description

4 MHz waterproof probe - vascular

For Sonotrax doppler (codes 29490 - 29491 - 29493).