



Product Code	29522
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279295221
PARAF	903141467

## Description

### FOETAL MONITOR PAPER

- Manufacturer of foetal monitor: BIONET
- Model: Fetalcare FC/700
- Standard FHR rate/BPM: US 60/240
- Type: roll
- Width: 215 mm
- Length: 23 m
- Grid: green
- Minimum order: 5