


SINGLE PATIENT FOETAL TRANSDUCER BELTS 6x120 cm - pink-light blue



Product Code	29544
--------------	-------

Unit of sale	100 kits of 2 pcs
--------------	----------------------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279295443
---------	---------------

Description

SINGLE PATIENT FOETAL BELTS - 100 pairs (1 pink +1 light blue)

Single patient elastic foetal belts to avoid cross-contamination.
Soft and comfortable belts resistant to abrasion, corrosion, light and heat.
Great elasticity and insulation.

Performance: more than 1,000 stretches.

Easy to wash, made of 85% polyester and spandex.

Supplied with buttons.

Size: 120x6 cm

Thickness: 0.8 mm \pm 0.1 mm

Equally spaced buttonholes to adjust belt according to the needs.

COMPATIBILITY: Bionet, EDAN, Philips M1350, Philips FM 20 series, GE, Corometrics, Sonicaid, Huntleigh and most foetal monitors.