

REUSABLE FOETAL TRANSDUCER BELTS 6x150 cm - grey



Product Code 29545

Unit of sale box of 10 pcs

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279295450

Description

REUSABLE FOETAL TRANSDUCER BELTS - box of 10

Reusable elastic fetal belt resistant to abrasion, corrosion, light and heat.
Soft, comfortable and easy to wash and dry.
Performance: more than 100,000 stretches.

Made of 85% polyester and spandex.

Supplied with buttons.

Size: 150x6 cm

Thickness: 1.0 mm \pm 0.1 mm

Equally spaced buttonholes to adjust belt according to the needs.

COMPATIBILITY: Bionet, EDAN, Philips M1350, Philips FM 20 series, GE, Corometrics, Sonicaid, Huntleigh and most foetal monitors.