

FOETAL MONITOR PAPER 143x150 mm - pack - white paper



Product Code	29580
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2322675
CND	Z1302
EAN/UPC	8023279295801

Description

FOETAL MONITOR PAPER

- Manufacturer of foetal monitor: SONICAID
- Model: Meridian/800 - Team
- Standard FHR rate/BPM: White paper
- Type: pack
- Width: 143 mm
- Length: 150 mm
- Sheets: 300
- Grid: white
- Minimum order: 10