



Product Code	29600
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2415998
CND	Z12020703
EAN/UPC	8023279296006
PARAF	903141695

Description

COLPY - GIMA LED COLPOSCOPE

Easy to use LED colposcope mounted on a 5 wheels base (Ø 41 cm).

Made in Italy.

Technical Specifications

Magnifying: 8.2X with standard 10X eyepieces
 Field: 35 mm
 Optic system: stereo with inter-pupil adjustment and optic head at 45°
 Illumination: LED bulb
 Working distance: 28 cm
 Interpupillary distance: 48-77 mm
 Green filter
 Operating voltage: 100-240V, 50-60 Hz
 Norms: CLASS I Type B - IEC 601-1 CE 93/42/EEC.