



Product Code	29603
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2416224
CND	U0899
EAN/UPC	8023279296037

Description

DL1 DIGITAL - USB 2.0 CAMERA with grabbing software for Colpy

Fitted on one eyepiece, allows to show image on PC through USB connection, becoming a high definition videocolposcope. Includes camera connector.

Technical Specifications

CMOS 1/3" sensor
 Max resolution pictures 1280x1024 pixels
 Max resolution movie 440x480 pixels