



Product Code	29612
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1709597
CND	Z12020703
EAN/UPC	8023279296129

Description

ALLTION LED COLPOSCOPE - 9X

High quality and complete colposcopes, with excellent optics and cold, longlife white light. Compact, lightweight and easy to use offers superior quality at a competitive price. A special universal Teflon® ball joint gives user 360° easy motion for gross focus. Comes assembled, ready for use.

Manual: GB, FR, IT, ES.

Technical Specifications

Magnification 9x
 Focus length: 316
 Field of view: Ø >22 mm
 Depth of field: >8 mm
 Light intensity: > 25,000 lux
 Interchangeable adjustable eyepiece 12.5x
 Illumination: longlife 3.3V / 10W LED
 Superior optics with individually adjustable eyepiece
 Built-in rheostat for brightness adjustment
 Easily changeable, LED bulb
 User-selectable built-in green filter
 Two (gross/fine) focus adjustment handles
 Two (gross/fine) height adjustment handles (838-1194 mm)
 Smooth action control
 Power supply: A C 100-240 V - 50/60 Hz; DC 12V 3A