

ALLTION LED COLPOSCOPE - 3.75X, 7X, 15X > 28,000 Lux



Product Code	29613
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1709597
CND	Z12020703
EAN/UPC	8023279296136

Description

ALLTION LED COLpoSCOPE 3,75x, 7x, 15x - > 25,000 LUX

High quality and complete colposcopes, with excellent optics and cold, longlife white light. Compact, lightweight and easy to use offers superior quality at a competitive price. A special universal Teflon® ball joint gives user 360° easy motion for gross focus. Comes assembled, ready for use.

Manual: GB, FR, IT, ES.

Technical Specifications

Magnification 3,75x, 7x, 15x
Focus length: 300 mm
Field of view: 9.5/18.5/38 mm
Depth of field: 5/8/34 mm
Light intensity: > 25,000 lux
Interchangeable adjustable eyepiece 12.5x
Illumination: longlife 3.3V / 10W LED
Superior optics with individually adjustable eyepiece
Built-in rheostat for brightness adjustment
Easily changeable, LED bulb
User-selectable built-in green filter
Two (gross/fine) focus adjustment handles
Two (gross/fine) height adjustment handles (838-1194 mm)
Smooth action control
Power supply: A C 100-240 V - 50/60 Hz; DC 12V 3A