

AMNIOCENTESIS NEEDLE 21 G x 150 mm - sterile



| | |
|---------------|----------------|
| Product Code | 29641 |
| Unit of sale | box of 25 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|---------------|
| RDM (NSIS) | 21094 |
| CND | A010204 |
| EAN/UPC | 8023279296419 |

Description

AMNIOCENTESIS NEEDLES - STERILE

Needle 21 G x 150 mm

For transabdominal sampling of amniotic fluid (usually between 15th and 18th pregnancy week) for prenatal diagnosis of genetic disorders of the foetus. It is also used for cytogenetic and biochemical studies. Accurate sharpening of the needle allows an easy and atraumatic insertion. Its mandrel is provided with an ultrasound detectable tip allowing an echo-guided insertion, and it is colour coded for an easy identification of the needle size. The needle is provided with an adjustable stopper to determine insertion depth, and with a female luer lock connector to collect the amniotic fluid through a syringe.

Made in Italy.