

TRUGLYDE ABSORB. SUTURE gauge 3/0 circle 1/2
needle 22mm - 75cm - undyed



Product Code	22967
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	1781300
CND	H010101020103
EAN/UPC	8023279229677
GMDN	13908

Description

TRUGLYDE - braided & coated polyglycolic acid

- Gauge: 3/0
- Metric: 2.0
- Colour: Undyed
- Circle: 1/2
- Needle profile: round body
- Needle length: 22 mm
- Suture length: 75 cm

Technical details:

- braided and coated polyglycolic acid
- coating of calcium stearate and polycaprolactone
- sterilized by ethylene oxide
- shelf life of 5 years
- wound Support : 28-30 days
- mass Absorption: 70-90 days

Applications:

- for soft tissue approximation during general, gynecological and orthopedic surgeries

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.