



Product Code	29737
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2852638
CND	U0899
EAN/UPC	8023279297379
PARAF	977813056
GMDN	42537

Description

Sterile GIMA collector

- Material: plastic
- Shaft length: 19.5 cm
- Packed in a single sterile pouch in box of 100 pouches

With Wallach papette all it takes is a simple rotating motion to collect a total sampling endocervical and ectocervical cells, including the transformation area, with less trauma.