

SMALL COTTON SWABS - cotton 5 mm - sterile



Product Code	29753
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2987343
CND	U0899
EAN/UPC	8023279297539
GMDN	33722

Description

Sterile Cotton Swab

- Sterile
- Shaft material: plastic
- Ø Cotton: 0.5 cm
- Shaft length: 15 cm

Packed in box of 100 pcs.