

# SMALL COTTON SWABS - cotton 5 mm (10 plastic bags of 100)



Product Code	29754
Unit of sale	box of 1000 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2733131
CND	M9099
EAN/UPC	8023279297546
PARAF	903142038
GMDN	33722

## Description

### Cotton Swab

- Shaft material: wood
- Ø Cotton: 0.5 cm
- Shaft length: 15 cm

Packed in plastic bags of 100 pcs., in carton of 1000 pcs.