

# S.S. CUSCO SPECULUM-MEDIUM insulated for HF applications



|               |                |
|---------------|----------------|
| Product Code  | 29763          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2913960        |
| CND        | U0899          |
| EAN/UPC    | 8023279297638  |
| PARAF      | 903142089      |
| GTIN       | 08023279297638 |
| GMDN       | 35352          |

## Description

HF Cusco stainless steel speculum - medium

With insulation for HF applications and suction connector.