

PIPELLA 1 HOLE - sterile



Product Code	29826
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2853730
CND	U0899
EAN/UPC	8023279298260

Description

Pipella 1 hole - disposable, sterile - box of 100

For endometrial cytologic sampling.
It allows a painless and easy routine screening of the endometrium avoiding dilation and traction of the uterus neck.

Hole: Ø 2.0-2.5 mm
External diameter: 3.1 mm
Length of cannula: 23.5 cm
Length of handle: 26.5 cm

Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.