

# GYN&CUBE PESSARY diam. 28/38 mm, h 38 mm - standard



Product Code	29887
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

EAN/UPC	3662914003453
PARAF	977813094

## Description

Gyn&Cube pessary  $\varnothing$  28/38 mm, h 38 mm - standard

Silicone, autoclavable, latex free pessary with pyramidal shape and draining holes.  
Kevlar reinforced lanyard.

Multilanguage instructions: GB, FR, IT, ES, DE.

### Indications:

- stress urinary incontinence during sustained exertion (aerobics, tennis, jogging)
- moderate or exteriorized prolapse
- moderate cystocele and rectocele
- after childbirth, transient prolapse or incontinence