

GYN&CUBE PESSARY diam. 32/44 mm, h 44 mm - large



Product Code	29888
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

EAN/UPC	3401564583782
PARAF	977813106

Description

Gyn&Cube pessary ø 32/44 mm, h 44 mm - large

Silicone, autoclavable, latex free pessary with pyramidal shape and draining holes.
Kevlar reinforced lanyard.
Multilanguage instructions: GB, FR, IT, ES, DE.

Indications:

- stress urinary incontinence during sustained exertion (aerobics, tennis, jogging)
- moderate or exteriorized prolapse
- moderate cystocele and rectocele
- after childbirth, transient prolapse or incontinence