

ETHICON ETHILON MONOFILAMENT SUTURES - gauge 8/0 needle 6.5 mm



Product Code	22299
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2561872
CND	H010201010101
EAN/UPC	8023279222999
GMDN	38000

Description

ETHILON - non absorbable sterile monofilament sutures

- Ethicon code: 2808G
- Gauge: 8/0
- Metric: 0.4
- Colour: black
- Circle: 3/8
- Needle lenght: 6.5 mm
- Needle: BV130-5
- Suture lenght: 13 cm
- Minimum order: box of 12 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.