

# CUPS HYSTEROGR. 35 mm. - sterile



|               |                |
|---------------|----------------|
| Product Code  | 29922          |
| Unit of sale  | 1 pc           |
| Minimum order | 10             |
| Type          | Medical device |
| Class         | I S            |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2590118       |
| CND        | U080102       |
| EAN/UPC    | 8023279299229 |
| PARAF      | 903143182     |

## Description

### CUPS FOR HYSTEROGRAPHY - 35 mm - sterile

Rigid, transparent tapered cup with a central conical cannula. Equipped with a 25 cm long flexible, pliable manipulator joined to the cup and two tubes equipped with clamps. One of these provides the vacuum, while the other injects the radiopaque medium. Perfect adhesion to the cervix by selfcentering is ensured by the tapered shape of the cup and the progressive application of a vacuum to the cup. The suction cup/cannula unit is radiotransparent and allows the proximal area of the endocervix to be controlled. The radiologist is able to manipulate at a distance from the radioactive field while at the same time comfortably immobilizing the uterus. Latex free.