



Product Code	29940
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	No

RDM (NSIS)	1760184
CND	U080101
EAN/UPC	8023279299403
GMDN	32655

Description

ENDOGYN - DISPOSABLE, STERILE - BOX OF 25

For cytologic endometrial sampling. It allows a painless and easy routine screening of the endometrium avoiding dilation and traction of the uterus neck. An outer sheath covers a mandrel with asymmetrical "V" shaped blades on the top. Each arm features fine toothings and thin cutting edges which are able to remove and keep more endometrial tissue. The head of one of the two arms features a hemispherical head which, when it goes back inside the sheath, allows for previously removed tissue to be aspirated.

Asimmetrical V shape blades tip with protrusions for sampling

External Ø 3.0 mm - Length 20 cm.

Hysterometry can be measured by means of the reference points located at 3, 6 and 9 cm along the sheath

Endometrial sampler for cytologic examination