

# UMBILICAL CORD CLAMP - sterile



Product Code	29941
Unit of sale	box of 1000 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	Yes

RDM (NSIS)	3000458
CND	V0202
EAN/UPC	8023279299410

## Description

UMBILICAL CORD CLAMP - sterile - box of 1000 (10 boxes of 100)

Sterile umbilical cord clamp is used to securely close off the umbilical cord of newborns immediately after birth.

Single use.

Made of PE.

Individually packed.

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.