

# X-SAFE UMBILICAL CORD CLAMP - sterile



Product Code	29942
Unit of sale	box of 500 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	1382691
CND	V0202
EAN/UPC	0373011105003
GMDN	32192

## Description

### X-SAFE UMBILICAL CORD CLAMP - sterile

Disposable umbilical cord clamp to allow a proper clamping of all umbilical cords and prevent their outcoming from the closing comb. Individually packed in coupled cellophane (nylon/polyethylene) and medical paper to ensure safe hygienic conditions to user. Sterilized by ethylene oxide.

Extra safe clamp avoid accidental re-opening. Can be re-opened only by cutting back side.

Made in Europe.