

LUBRICANT GEL - tube 42 g - sterile



Product Code	29958
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1707200
CND	V9099
EAN/UPC	8698712459463
PARAF	976329639
GMDN	60796

Description

KLY LUBRICATING GEL

- One 42 g tube
- Sterile and paraben-free

Hypoallergenic lubricating gel, does not irritate skin, dissolves in water and is easily cleaned.

It does not contain oil or oily substances.
It has no toxic effects or odour.

It is used in endoscopy, gastroscopy, colonoscopy and rectoscopy.

Multilanguage box and manual: GB, FR, IT, ES, DE, PT, GR, RU, TR, Arabic.