

CUSCO SPECULUM STERILE - medium



|               |                |
|---------------|----------------|
| Product Code  | 29977          |
| Unit of sale  | box of 120 pcs |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I S            |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |                |
|------------|----------------|
| RDM (NSIS) | 414971         |
| CND        | U089006        |
| EAN/UPC    | 8023279299779  |
| GTIN       | 08023279299779 |
| GMDN       | 37468          |

Description

CUSCO Speculum - medium 90 x 28 mm

Disposable speculum with gradual screw-action adjustment. It is manufactured with the best non-toxic transparent material, comes in a single transparent vacuum pack and is sterilized with ethylene oxide.

Sterilization is guaranteed by the special quality control label on each pack.

Available 3 sizes distinguishable through screw colour.  
Made in Europe.