


SILICONE RING PESSARY WITH SUPPORT diameter 102 mm



Product Code	30049
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2770487
CND	U089005
EAN/UPC	6972994800208

Description

Silicone Ring Pessary with Support - white Ø 102 mm

Used for first and second degree prolapses that are complicated by a mild cystocele. Single patient. Individually packed in a cardboard box.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, SE, LV, PL, HU, SI, HR, RO, BG, GR, Arabic.

Pessary is used to support uterine prolapse in women, including cystocele and rectocele, as well as stress urinary incontinence. This medical device provides instant relief, allowing the patient to resume her normal daily activities immediately. The device must be prescribed and positioned by a qualified specialised practitioner.

Latex free. Silicone material ensures softness and flexibility.