

# SILICONE GELLHORN PESSARY WITH LONG STEM

## diameter 64 mm



|               |                |
|---------------|----------------|
| Product Code  | 30054          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2770490       |
| CND        | U089005       |
| EAN/UPC    | 6972994800253 |

### Description

Silicone Gellhorn Pessary with Long Stem - white Ø 64 mm

Used for a second to third degree uterine prolapse or procidentia. Single patient. Individually packed in a cardboard box.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, SE, LV, PL, HU, SI, HR, RO, BG, GR, Arabic.

Pessary is used to support uterine prolapse in women, including cystocele and rectocele, as well as stress urinary incontinence. This medical device provides instant relief, allowing the patient to resume her normal daily activities immediately. The device must be prescribed and positioned by a qualified specialised practitioner.

Latex free. Silicone material ensures softness and flexibility.